

Business from technology

Pneumatically powered material testing devices for the extreme test conditions

Pekka Moilanen, Timo Saario and Jussi Solin Advancing Materials testing in Hydrogen Gas Sandia National Laboratories April 9-10, 2013





Our strategy: Materials performance in simulated extreme conditions, real failure mechanism with real stressors in real environment.

VTT Materials research



VTT Materials Performance Research Hall 1 in Espoo (altitude: -40m)




- Background
- The pneumatic servo controlled material testing system
 - Applications and co-operators
- In-situ tensile and fatigue tests in SCK-CEN
 - In-situ dynamic tests vs. post irradiation tests
- In-situ fuel cladding device, MeLoDie, CEA
 - Stress, strain variations during the irradiation
- Miniature size of the autoclave for SCW environment
 - Double bellows loading apparatus, tests at LWR and SCW
- Liquid lead testing device in JRC Petten
 - Conceptual design for LL testing devices
 - On-line load monitoring system with the pneumatic double² bellows load apparatus
 - Load line displacement measurement system for crack growth measurement

→ Hybello Hydrogen testing system in VTT



Pneumatic system;


- No moving parts on the pressure boundary
 → more sensitivity, the minimal leakage.
- Possible to integrate into a very small test chamber.
- Possible to generate the test load with an inert gas (Helium).





The pneumatic loading unit;

- Test load generation with axial movement
- $F=p_{in}/p_{out}*A_{eff}$
- p_{in}= internal pressure
- p_{out}= external pressure
- A_{eff}= effective cross section (calibration)
- Axial movement ±2 mm
- The maximum load 50 kN
- The minimum load 100 N




Background



The load frame


- -Support for the specimen, LVDT sensor and pneumatic loading unit
- -Alignment for the specimen
- -Depending on test type, specimen type and size
- -In-reactor creep fatigue device






Control setup

- -Programmable servo controller






Background

a.sji.enable 1 a.sji.inptrO ai2.outvalue a.fi.inptr a.sji.outvalue a01.inptr a.fi.outvalue a.sji.inptr1 upar1 timer1.enable 1 fnkey 4 timer fnkey 6 timer conflag 1 runflag 1 relay 1

The air booster based Helium pressure adjusting loops

Background




Designed for use in OSIRIS-reactor in France



Development of the calibration systems for the pneumatic loading unit


Background




Background

1. Determination of the own stiffness of the bellows

OS = 3.15 bar/mm





2. Determination of the effective cross section for the bellows

$$A_{eff} = 86.5 \text{ mm}^2$$



In-reactor experiment 2000-2007







The pneumatic tensile test modules which were designed and constructed at VTT have been installed in the BR-2 reactor at SCK-CEN, Belgium.

First unique experiments i.e., strain-controlled tensile experiments using OFHC-Cu have been carried out in-situ under neutron irradiation.

Irradiation conditions:

BR-2 position G60 reactor pool water at 90° C neutron flux $\sim 0.3 \times 10^{14} \text{ n cm}^{-2}\text{s}^{-1}$ (E>1 MeV) damage rate $\sim 2 \times 10^{-4} \text{ dpa h}^{-1}$


Loading conditions:

uniaxial tensile test strain controlled strain rate ~10⁻⁷ s⁻¹

European Fusion Technology Programme: VTT, Risø, SCK-CEN




The dynamic response of a tensile specimens during un-irradiated, post-irradiated and in-reactor experiences







In-core creep fatigue testing device



Tests carried out at Mol test reactor BR-2 in 2003-2005



In-core multi-axial creep test device

- VTT CEA in-kind project
- Bi-axial loading; push and pull bellows
- Tubular pressurized specimen
- Constant load tests with axial strain/hoop strain sweeps
- He-re-circulation loop



Osiris test reactor (CEA)





First irradiation (OSIRIS-reactor, France) in 2013



Miniature size crack growth rate testing device


- Developed together with VTT, Finland and JRC Petten, The Netherlands (2005...on going)
- Supercritical water (SCW) environment: 650 °C and 350 bar
- Autoclave outer diameter: 64 mm
- Autoclave inner diameter: 32 mm
- The maximum load: ~±3 kN
- Displacement range ± 1 mm
- Test types: Constant load, rising load, fatigue






D:\Project\VTT\Autoklaavi\FEM\Autoklaavi_1-3.m

RESULTS: 4- B.C. 2.DISPLACEMENT_1.OPERATION DISPLACEMENT - MAG MIN: 3.99E-06 MAX: 2.32E+00 DEFORMATION: 4- B.C. 2.DISPLACEMENT_1.OPERATION DISPLACEMENT - MAG MIN: 3.99E-06 MAX: 2.32E+00 FRAME OF REF: LOCAL 6.






30/04/2013



Miniature size crack growth rate testing device

- SCW environment
- Double bellows loading apparatus
- Test set-up for the crack growth test with 5 mm DC(T) specimen under SCW coolant conditions





VTT-HyBello: pneumatically powered test device for hydrogen gas, ready in 2013

Target test:

Specimen: 7≤Ø≤8mm Drilled: 100≤Ø≤500µm

Nf criterium: Crack~500µm

Instrumentation: $\Delta a/\Delta N$

Target Equipment:

H₂: ≤350bar,≤100°C

Axial loading: ±30kN (~1%) Internal load train: 2 bellows Frictionless load: Δp (p123)


Medium: air, He or water Frequency: LCF -> HCF

Safety concern:

Vessel volume: ≤0.5 liter H₂ volume: << 0.5 liter

Infra: 220V, gas supply -> find safe

place/H2 leaks, water pool?





VTT-HyBello: pneumatically powered test device for hydrogen gas, ready in 2013

adapter


Strain amplifier

Target in EU project MATHRYCE:

Specimen: 7≤Ø≤8mm


Drilled notch: 100≤Ø≤500μm N_f criterium: crack a~500μm

To measure: $\Delta a/\Delta N$



Some day ? Strain controlled LCF tests:

- Cyclic stress strain response (CSS)
- ASME III type desing curve








FABELLO for ϕ 8 mm DIRECT STRAIN CONTROL

(Results in) (PVP 2013) → ASTM E 606 valid LCF data in 135 bar 325 °C water



This in 1000 bar Hydrogen ≈ ??




Safety concerns and open questions

- Amount of H₂ gas in the test section:
 Possible solution: the minimum volume of the H₂ vessel and small size pressure tubes; Dout = 1.6 mm Din = 0.7 mm
- Leakage of H₂:
 Possible solution: Water pool, safe test environment +
 Leagage detectors -> ventilation of the test room
- The strain measurement system (electrical power is needed for the LVDT sensor or strain gauges) -> sparking?
 Possible solution: place the LVDT sensor outside from the H₂ test section
- Sealing of the autoclave and needed feedthroughs:
- Solution: Double coinical metal rign sealing, Conax or Swagelok type of the standard connectors for the wires.
- Material selection (alloy 625), coating, welding etc...

30/04/2013



Illustration for the H2 test set-up





High precision Pneumatic material testing technology

- No moving parts on pressure boundary
- Patented pneumatic servo controlled pressure adjusting loop
 - ⇒ more sensitivity, more flexibility
- Easy to move and integrate into different test environments
 - Autoclave testing at laboratory
 - Tensile/fatigue tests inside the test reactor core with 35 m pipe lines
- Possible to test many specimens at the same time
- Possible to test with inert gas (Helium) ⇒ no pollution problems for the test environment
- Pneumatic double²bellows load apparatus -> direct load measurement under high temperature and pressurized test environments